

IHE Integration Statement			
Vendor	Product Name	Version	Date
Corepoint Health	Corepoint Integration Engine	7.2	2/18/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Audit Trail and Node Authentication	Secure Node	None
Consistent Time	Time Client	None
Healthcare Provider Directory	Provider Information Consumer	None
Mobile Alert Communication Management	Alert Aggregator	None
Mobile Care Services Discovery	Care Services Selective Consumer	None
Mobile Care Services Discovery	Care Services Selective Supplier	None
Mobile Care Services Discovery	Care Services Update Consumer	None
Mobile Care Services Discovery	Care Services Update Supplier	None
Mobile access to Health Documents	Document Consumer	None
Mobile access to Health Documents	Document Recipient	None
Mobile access to Health Documents	Document Source	None
Patient Demographics Query	Patient Demographics Consumer	None
Patient Demographics Query	Patient Demographics Supplier	None
Patient Demographics Query for Mobile	Patient Demographics Consumer	None

Internet address for vendor's IHE information :
<http://www.corepointhealth.com/sites/default/files/integration-statement-corepoint-integration-engine.pdf>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Corepoint Health	Corepoint Integration Engine	7.2	2/18/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Patient Demographics Query for Mobile	Patient Demographics Supplier	None
Patient Demographic Query HL7 V3	Patient Demographics Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-referencing for MPI	Patient Identity Cross-reference Manager	None
Patient Identifier Cross-reference for Mobile	Patient Identifier Cross-reference Consumer	None
Patient Identifier Cross-reference for Mobile	Patient Identity Cross-reference Manager	None
Patient Identifier Cross-Reference HL7 V3	Patient Identifier Cross-reference Consumer	None
Query for Existing Data for Mobile	Clinical Data Consumer	None
Restricted Metadata Update	Update Initiator	None
Cross-Community Access	Initiating Gateway	None
Cross-Community Access	Responding Gateway	None
Cross Community Patient Discovery	Initiating Gateway	None

Internet address for vendor's IHE information :
<http://www.corepointhealth.com/sites/default/files/integration-statement-corepoint-integration-engine.pdf>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net

IHE Integration Statement			
Vendor	Product Name	Version	Date
Corepoint Health	Corepoint Integration Engine	7.2	2/18/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Cross Community Patient Discovery	Responding Gateway	None
Cross-Enterprise Document Reliable Interchange	Document Source	None
Cross-Enterprise Document Sharing	Document Source	None

Internet address for vendor's IHE information :
<http://www.corepointhealth.com/sites/default/files/integration-statement-corepoint-integration-engine.pdf>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net